


 **Ionizer Replacement Instructions**

PARTS LIST:

- Acrylic Top Plate (1/4" x 170 mm x 170 mm) (1x)
- 100078 Makerbeam 10 mm x 10 mm x 100 mm (4x)
- 100179 Makerbeam 10 mm x 10 mm x 150 mm (4x)
- 100988 Makerbeam Corner Cube (4x)
- 92125A126 18-8 Stainless Steel Hex Drive Flat Head Screw, M3 x 0.5 mm Thread, 6 mm Long (8x)
- 92125A136 18-8 Stainless Steel Hex Drive Flat Head Screw, M3 x 0.5 mm Thread, 20 mm Long (4x)

MATERIAL 1. Aluminum 2. Acrylic	UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES.	 Stanford Research Systems		
	TOLERANCES: (UNLESS OTHERWISE NOTED) .X ± .XX ± .05 .XXX ± ANGLE ±	1290-D Reamwood Ave, Sunnyvale CA 94089		
FINISH 2.		TITLE RGA Work Stand		
APPROVED	SIZE	CODE	SRS PART NO.	REV
CHECKED	B			
DRAWN A.B.	12/20/2024	SCALE 1:2	WEIGHT	SHEET 1/1