


REVISIONS			
REV	DESCRIPTION	DATE	APPROVED
B	CHANGED INSIDE HOLE DIAMETER CHANGED TOLERANCE.	4/25/22	GAW



DETAIL A
SCALE 5.000
DIMENSIONS ARE TYPICAL

DETAIL B
SCALE 4.000
3 EQUAL SPACES

- NOTES: UNLESS OTHERWISE SPECIFIED
- MATERIAL: STAINLES STEEL COIL 302/304 HALF HARD.
.010 THICK.
 - FINISH: BLACK OXIDE.
 - REMOVE ALL BURRS AND SHARP EDGES.
 - PART MARKING: BAG AND TAG WITH VENDOR ID,
MONTH/YEAR OF MANUFACTURE PART NUMBER AND
CURRENT REVISION LEVEL.

MATERIAL 1	UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES. TOLERANCES: .X ± .XX ± .XXX ± .005 ANGLE ±	STANFORD RESEARCH SYSTEMS, INC. 1290-D Redwood Avenue, Sunnyvale CA 94089	
FINISH 2		DRAWING TITLE: CHOPPER BLADE, VARIABLE DUTY CYCLE 10-90 6 SECTOR, SR542	
ENG BY/DATE A BERGER	SIZE/PROV/E FILENAME: D	SRS PART NO. 7-02831-716	REV B
DWN BY/DATE G.WATERS	Jun-04-20	SCALE: 4.00	SHEET: 1 OF 1

NEXT ASSY	USED ON	DO NOT SCALE DRAWING
-----------	---------	----------------------