



Stanford Research Systems, Inc.
1290 C Reamwood Avenue, Sunnyvale, CA 94089, U.S.A.
Phone: 408-744-9040 FAX: 408-744-9049
Email: janied@thinkSRS.com .. Website: www.thinkSRS.com

4/18/19

Volatility Statement for SRS BGA244 (Binary Gas Analyzer)

The BGA244 contains two types of memory, as listed below:

Non-Volatile Memory:

256 Mbits Flash Contains instrument's firmware, as well as user modifiable instrument settings.

To clear all user settings:

Recall "Factory" setting from the Store|Recall menu.

Volatile Memory:

256 Mbits SDRAM Used for operation

To clear:

Cycle the power

Please feel free to contact SRS should you have any further questions.

A handwritten signature in black ink, appearing to read 'Janie Du', written in a cursive style.

Janie Du
Sales & Marketing