

Certificate of Volatility				
<b>Model:</b> SIM928	<b>Part Number:</b>	<b>Manufacturer:</b> Stanford Research Systems		
		<b>Street Address:</b> 1290-D Reamwood Ave		
		<b>City:</b> Sunnyvale	<b>State:</b> CA	<b>Zip:</b> 94089
<b>Item Nomenclature:</b>				
Volatile Memory				
Does the item contain volatile memory (i.e., memory whose contents are lost when power is removed)? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No				
If the answer is 'Yes', please provide the following information for each type (use additional sheets if required):				
Type (SRAM, DRAM, etc.): SRAM	Size: 1 kB	User Modifiable: <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	Function: Firmware general storage and user communications data buffers	Process to Clear: Remove power
Type (SRAM, DRAM, etc.):	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Type (SRAM, DRAM, etc.):	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Non-Volatile Memory				
Does the item contain non-volatile memory (i.e., memory whose contents are retained when power is removed)? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No				
If the answer is 'Yes', please provide the following information for each type (use additional sheets if required):				
Type (BBRAM, Flash, EEPROM, etc.): EEPROM	Size: 768 bytes	User Modifiable: <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	Function: Factory calibration data, and user configuration	Process to Clear: send remote "*"RST" command before removing power. Factory cal data cannot be modified by user
Type (BBRAM, Flash, EEPROM, etc.): Serial EEPROM	Size: 2 kbyte	User Modifiable: <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	Function: records battery pack usage	Process to Clear: remove battery pack
Type (BBRAM, Flash, EEPROM, etc.): Flash	Size: 32 kbyte	User Modifiable: <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	Function: program image	Process to Clear: remove CPU -- destructive to instrument
Media				
Does the item contain media storage capability (i.e., removable or nonremovable disk drives, tape drives, memory cards, etc.)? <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No				
If the answer is 'Yes', please provide the following information for each type (use additional sheets if required):				
Type (Disk, Tape, etc.):  Removable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Type (Disk, Tape, etc.):  Removable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Type (Disk, Tape, etc.):  Removable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
<b>Additional Information:</b>				
<b>Vendor/Manufacturer Representative Information</b>				
<b>Name:</b> Janie Du	<b>Title:</b> Sales	<b>Office Phone:</b> 408-744-9040	<b>Fax/Email:</b> 408-744-9049/ info@thinkSRS.com	