



Stanford Research Systems, Inc.

1290 C Reamwood Avenue, Sunnyvale, CA 94089, U.S.A.

Phone: 408-744-9040 FAX: 408-744-9049

Email: janied@thinkSRS.com .. Website: www.thinkSRS.com

March 21, 2019

Volatility Statement for SR630

The SR630 contains a Static RAM which holds calibration bytes and instrument settings. It is backed up by a 3V lithium battery.

Location	Type	Size	Function	User Modifiable?	Data Input Method	To clear User Settings
U 102	Static RAM	32K x 8 70 nS	Instrument Settings and Cal Bytes	Yes	Front Panel RS-232 GPIB	Power ON with BKSP * Or, remove the battery (BT100)

* With power off, holding down BKSP button while powering on the instrument (clear instrument settings and bring back factory default cal bytes.)

Please feel free to contact SRS should you have any further questions.

Janie Du
Sales & Marketing